



DEVICE PARTS	INTENDED PURPOSE	
	Microperforator for aspiration of drugs to be used with crn®CRONO®Reservoir syringes.	
	CLASSIFICATION	
	Class Is medical device, according to Rule 2 of Annex IX of 93/42/EEC European directive and subsequent changes introduced by the 2007/47/EC European directive.	
	MANUFACTURER	
	MULTIMEDICAL s.r.l. Via Guido Rossa, 71 - Zona Industriale Gerbolina 46019 Viadana (MN) - Italy Tel.: +39-0375-785882	
	CE MARKING	 0123
	ITALIAN REPERTORY NUMBER OF MEDICAL DEVICES (RDM)	1364422
	ITALIAN CLASSIFICATION OF DEVICE (CND)	A0704
	CANÈ REF	SPIKE/04
MANUFACTURER REF	091548/1	

PART No.	DESCRIPTION	MATERIAL
1	Microperforator with female Luer Lock connector (ISO 80369-7)	ABS
2	Tip cover	LDPE
3	Air filter	PP 3 µm hydrophobic membrane
4	Male Luer-Lock cap (ISO 80369-7)	ABS

DEVICE SPECIFICATION
Sterile, single-use, non-pyrogenic, non-resterilizable product. Device free of latex and phthalate components (DEHP).

PACKAGING	GTIN	PCS	FEATURES	MATERIAL
Primary	08050046231373	1	Blister 12x5 cm	Medical degree Paper PP/PE Plastic film
Secondary	18050046231370	60*	Bag	PE
Tertiary	28050046231377	300	Box 39x29x20 cm	Corrugated cardboard

* each bag of 60 pieces contains an instruction sheet for use.

ENVIRONMENTAL CONDITIONS FOR STORAGE AND TRANSPORT	EXPIRY DATE
Store the product away from direct exposure to sunlight, in dry and clean environments and away from extreme temperatures and excessive humidity.	4 years and 11 months from the date of manufacture.

 **WARNINGS**

Disposable consumables must be disposed of in accordance with the regulations governing the disposal of medical waste.