


DEVICE PARTS	DESCRIPTION	
	Disposable filtering device, usable with any type of infusion set.	
	CLASSIFICATION	
	Class Is medical device, according to Rule 2 of Annex IX of 93/42/EEC European directive and subsequent changes introduced by the 2007/47/EC European directive.	
	DISTRIBUTOR	
	CANÈ S.p.A. Via Cuorgnè 42/a 10098 Rivoli (TO) - Italy Tel.: +39-011-9574872	
	MANUFACTURER	
	MULTIMEDICAL S.r.l. Via Guido Rossa, 71 - Zona Industriale Gerbolina 46019 Viadana (MN) - Italy Tel.: +39-0375-785882	
	CE MARKING	CE 0123
	ITALIAN DEVICE REPERTORY REGISTRATION	88883/R
	ITALIAN DEVICE CLASSIFICATION (CND)	A03020101
REF CANÈ	FILTRAJET	
REF MANUFACTURER	090500	

PART N.	DESCRIPTION	MATERIAL
1	Luer Lock Male cap	LDPE
2	Female luer lock connection	ABS
3	Pediatric filter 0,2µ	Polysulfone membrane
4	Male luer lock connection	ABS
5	Labyrinthine Female Luer Lock cap	PE
6	2 transparent tubes	PU 93shA

DEVICE SPECIFICATION

Total length 17,5cm ± 0,5cm.
 Tube Ø: 0.5mm internal – 2.3mm external.
 Sterile, single-use, non-pyrogenic, non-resterilizable product.
 Device free of latex and PVC components.

PACKAGING	GTIN	PCS	FEATURES	MATERIAL
Primary	08050046230352	1	Blister cm 5x24	Medical paper PP/PE Film
Secondary	18050046230359	50	Cardboard box cm 28,5x15,8x7,5	Corrugated cardboard
Tertiary	28050046230356	500	Cardboard box cm 59x39x18	Corrugated cardboard

ENVIRONMENTAL CONDITIONS FOR STORAGE AND TRANSPORT
EXPIRY DATE

Store the product away from direct exposure to sunlight, in dry and clean environments and away from extreme temperatures and excessive humidity.

4 years and 11 months from the date of manufacture.


WARNINGS

Disposable consumables must be disposed of in accordance with the regulations governing the disposal of medical waste.