


DEVICE PARTS	DESCRIPTION
	<p>Device intended for the connection, via a tubular apparatus equipped with appropriate fittings, of medical devices, for the administration of infusion solutions/drugs/anesthetics or gases to the patient.</p>
	<p><b>CLASSIFICATION</b></p> <p>Class IIa medical device, CE marked according to Directive 93/42/EC and subsequent updates.</p>
	<p><b>MANUFACTURER</b></p> <p><b>MULTIMEDICAL s.r.l.</b> Via Guido Rossa, 71 - Zona Industriale Gerbolina 46019 Viadana (MN) - Italy Tel.: +39-0375-785882</p>
	<p><b>CE MARKING</b></p> <p>CE 0123</p>
	<p><b>ITALIAN REPERTORY NUMBER OF MEDICAL DEVICES (RDM)</b></p> <p>88981</p>
	<p><b>ITALIAN CLASSIFICATION OF DEVICE (CND)</b></p> <p>A03020101</p>
	<p><b>CANÈ REF</b></p> <p>Y-SET*MM</p>
	<p><b>MANUFACTURER REF</b></p> <p>090765</p>

PART No.	DESCRIPTION	MATERIAL
1	3 transparent tubes	PU 93shA
2	2 male Luer-Lock connectors	ABS
3	Female Luer-Lock connector	ABS
4	Y-Connector	ABS
5	2 Female Luer-Lock caps	PE
6	Male Luer-Lock cap	LDPE

DEVICE SPECIFICATION
<p>Total length 10 cm. Tube Ø: 0.5mm internal – 2.3mm external. Sterile, single-use, non-pyrogenic, non-resterilizable product.</p>

PACKAGING	GTIN	PCS	FEATURES	MATERIAL
Primary	08050046231465	1	Blister 4x16 cm	Medical degree Paper PP/PE Plastic film
Secondary	18050046231462	60*	Box 38x23x12cm	Corrugated cardboard

\* Each box of 60 pieces contains a sheet of instructions for use.

ENVIRONMENTAL CONDITIONS FOR STORAGE AND TRANSPORT	EXPIRY DATE
Store the product away from direct exposure to sunlight, in dry and clean environments and away from extreme temperatures and excessive humidity.	4 years and 11 months from the date of manufacture.

 **WARNINGS**

Disposable consumables must be disposed of in accordance with the regulations governing the disposal of medical waste.